

Merlin

Reading • Writing • Hobbies • Computer Compatible

The Clear Choice in Desktop Video Magnifiers



enhanced
vision[™]

Merlin

Whether you're at home, school or at the office, now you can regain your independence with Merlin. Read books, magazines and important documents. Write letters, view color photos, do crossword puzzles and even pursue a favorite hobby.

See for yourself what makes Merlin the greatest value and best picture quality available in desktop video magnification. Enjoy a feature rich system made to fit your needs and budget.



Easily pivot Merlin's LCD screen horizontally and vertically to provide a comfortable viewing position. View oversized documents and work on special projects with ease on Merlin's large workspace.

“ The Merlin is the best investment I have made for myself in a long time. Nobody can put a price tag on being able to read, write and see my mail again! ”


—Jeff Helms, CA



Plug in the new Feature Kit to add line markers, windowing/masking and VGA toggle to any Merlin LCD model.

Features Include:

- Simple, easy-to-use tactile controls
- Long lasting LED lighting for truer picture and color
- Ships in one box, simple to unpack, just plug in and begin
- Available in 19" & 22" monitors
- Screen easily pivots, both horizontally and vertically to provide the most comfortable viewing position
- Wide range of adjustable magnification from 2.4x to 77x (varies with display size)
- Computer Compatible (Plus Models or additional hardware required)
- 7 preset viewing modes
- Up to 28 custom programmable color select options available
- Slim, lightweight profile
- Adjustable brightness
- Auto focus
- 3 year warranty



Pivot and adjust the Merlin screen to maximize peripheral vision.

Merlin

Merlin's unique screen shifts from left to right and up and down for your most comfortable viewing position. Providing a clear crisp image and truer picture and color, you'll find reading and performing daily tasks enjoyable again.



Merlin Plus

Built for users who want a computer compatible system. Merlin's split screen mode lets you view magnified images and computer files simultaneously. Also enjoy horizontal and vertical line markers and windowing. Merlin Plus has it all.



To schedule a no obligation demonstration, please call (888) 811-3161.

Merlin VGA

Reading • Writing • Hobbies • Computer Compatible

Merlin VGA

Merlin VGA offers many of the same features of our standard Merlin in a low cost package, allowing you to easily connect your own VGA monitor for superior image quality and maximum versatility.

Merlin VGA Plus

Built for users who want a full-featured system, this computer compatible video magnifier has it all. View computer files and magnified images simultaneously.

Features Include:

- Affordable
- Connects to any computer monitor, LCD, or Plasma television with VGA inputs
- Intelligent Focus Lock Mode for writing
- Object Locator Mode
- 3 year warranty
- 7 preset viewing modes
- Up to 28 color select options available
- Computer Compatible (Merlin VGA Plus)
- Light On/Off Switch



Local Dealer:

enhanced
vision™

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

To schedule a no obligation demonstration, please call (888) 811-3161.