

amigo

High-Definition • Light • Portable

Your Desktop Magnifier On The Go


Amigo, your desktop video magnifier on the go, is light and simple to use. Place it directly on your reading material and adjust magnification with an easy turn of the dial. Read price tags, menus, labels and more, anywhere!

Amigo Features:

- Crisp, high-definition image
- Lightweight and Portable (only 1.3 lbs)
- 3.5x to 14x adjustable magnification
- Large 6.5" anti-glare LCD screen tilts for most comfortable viewing angle
- Large field of view allows the user to see more of the reading area
- Freeze frame with adjustable size and contrast
- Removable 2-hour battery with internal charger (spare battery and external charger also included)
- 6 viewing modes to optimize contrast and brightness
- Connects with any TV for increased magnification
- Writing stand and carrying case included
- 2 year warranty
- Also available in black



"The Amigo is worlds ahead of the typical handheld magnifier. It goes with me almost everywhere I go. I use it at all veteran related meetings..." - **Cornelius J. Appleby, PA**



enhanced
vision[™]

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

amigo

A portable magnifier for any task...



Write Checks



Read Menus



Read Price Tags



Read Ingredients

Local Dealer:

To schedule a no obligation demonstration, please call (888) 811-3161.

enhanced
vision[™]

(888) 811-3161
(714) 465-3400
www.enhancedvision.com