

pebble-mini

Ultra-Compact • Portable • Lightweight

The ultimate lightweight hand-held magnifier!

Pebble-mini's ultra-compact design makes it the ideal hand-held magnifier whether you're at home or on the go. Easily read prescriptions, mail, labels, price tags, menus and so much more.

Wear Pebble-mini around your neck or carry it in your pocket or purse – the only time you'll know it's there is when you need it!

Pebble-mini Features:

- Adjustable magnification from 2x to 10x
- Lightweight design (less than 4 ounces)
- 3.0" high resolution LCD display
- Large viewing area in a small package
- Adjustable brightness
- Large and easy-to-use controls
- 28 viewing modes available
- Freeze image feature with capability to magnify and change mode
- Image capture with save and recall
- Upload stored images to PC
- 2 to 3 hours of battery usage
- Charge with USB into PC or outlet
- Protective pouch and neck strap
- Hands-free reading stand
- 2 year warranty



“ I really enjoy how light the Pebble-mini is, it fits in my pocket and I can wear it around my neck when I'm shopping allowing me to see so much more! ”

– Nancy Becker



enhanced
vision

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

pebble-mini

Ultra-Compact • Portable • Lightweight



or



Standard
Optical
Magnifiers

Pebble-mini

- Adjustable Magnification
- 28 Viewing Modes Available
- Freeze Frame Feature
- Image Capture
- Fully Portable
- Large Field of View
- Adjustable Brightness



Local Dealer:

To schedule a no obligation demonstration, please call
(888) 811-3161

enhanced
vision
(888) 811-3161
(714) 465-3400
www.enhancedvision.com