

# transformer

USB Computer Compatible • Portable • Lightweight

## The next generation of portable electronic magnifiers:

Transformer is the most flexible and portable solution for reading, writing and viewing magnified images at any distance. Compatible with popular magnification software and the latest windows operating systems, Transformer's unique compact design is perfect for school, work, or home.

### Transformer Features:

- Weighs less than three pounds
- 2.4x to 30x magnification (based on 17" screen)
- Camera rotates 330° for reading, distance & self viewing modes
- Auto install software / plug & play
- USB connectivity for laptop/computer
- Compatible with popular magnification software & the latest windows operating systems
- Capture and store images
- 28 custom color select modes
- Self contained on/off LED lighting
- Battery operated - up to 4 hours
- Includes soft carrying case & cloth sleeve
- Memory setting for each camera position
- 2 year warranty



*The Transformer makes keeping up with classwork possible again!*



**enhanced**  
**vision**

**(888) 811-3161**

**(714) 465-3400**

**[www.enhancedvision.com](http://www.enhancedvision.com)**

# transformer



Transformer's sleek compact design and traveling case is ideal for easy transport.




To schedule a no obligation demonstration, please call (888) 811-3161.

Local Dealer:



Read, write, be independent again!



See near, far, and everywhere in-between

**enhanced**  
**vision**

**(888) 811-3161**

**(714) 465-3400**

**www.enhancedvision.com**